

EZ Clean Reamer Handle



The EZ Clean Reamer Handle allows for quick assembly and disassembly, simplifying the procedure as well as cleaning. The cross design ensures a stable connection with various reamer sizes.

A SIMPLE YET STRONG CONNECTION

- The cross design features a simple and stable connection of the reamer with no movement. This enables quick assembly and disassembly of the reamers, providing ease of use.
- EZ Clean disassembly is comprised of a simple “push and twist” of the locking collar, allowing all the parts of the cross connection to slide apart for cleaning.

RELIABILITY THROUGH PRODUCTIVITY

- Smaller head diameter creates better visibility with a smaller incision.
- A free floating Teflon sleeve over the shaft of the reamer handle protects soft tissues from becoming damaged.

Greatbatch Medical Headquarters
10000 Wehrle Drive
Clarence, NY 14031
tel 716-759-5600
fax 716-759-2562

US Customer Service
4532 Park 30 Drive
Columbia, City, IN 46725
tel 260-244-6300
fax 260-244-6345

Europe & Asia-Pacific Customer Service
L'Echelette 7
2534 Orvin
Switzerland
tel +41 32 358 01 11
fax +41 32 358 01 00

orthopaedicsna@greatbatchmedical.com

orthopaedicseurope@greatbatchmedical.com