

Handles



Greatbatch Medical offers a complete line of handle options. Whether you require silicone over-molding to a specific form, a precise torque-limiting handle or a ratchet handle, a wide selection is available.

DESIGNED WITH PERFORMANCE IN MIND

- Robust gear design and internal bearing provide smooth ratcheting performance.
- Ratchet handle thumb actuator allows for precise, one-handed actuation.
- Torque handle performance maintained within $\pm 10\%$ of nominal value over lifecycle.

RELIABLY BUILT TO YOUR NEEDS

- Solidly constructed of USP Class VI Silicone, over-molded on an aluminum or stainless steel core.
- Available design options fit standard industry shapes to allow for drop-in replacement.
- Choose a quick connect coupling or a custom shaft designed to your specifications.

Patent Pending

Greatbatch Medical Headquarters
10000 Wehrle Drive
Clarence, NY 14031
U.S.A.
tel +1 716-759-5600
fax +1 716-759-2562

US Customer Service
4532 Park 30 Drive
Columbia City, IN 46725
U.S.A.
tel +1 260-244-6300
fax +1 260-244-6345

Europe & Asia Customer Service
L'Echelette 7
2534 Orvin
Switzerland
tel +41 32 358 0111
fax +41 32 358 0100

orthopaedicsna@greatbatchmedical.com

orthopaedicseurope@greatbatchmedical.com

www.greatbatchmedical.com